Pantoprazole 40 mg Tablet

Structure:

Molecular Formula and Mass: C₁₆H₁₅F₂N₃O₄S - 383.37 g/mol

Category: Proton Pump Inhibitor

Sample:

Grind one tablet and dissolve in 120 mL of methanol. Shake at least 10 min and filter. 40.0 mg/120 mL = 0.333 mg/mL, which is the required concentration representing 100%.

Standards:

High Standard:

The high limit is 115%; therefore the concentration of the high standard = $(0.333 \text{ mg/mL} \times 1.15 = 0.383 \text{ mg/mL}$. Weigh approximately 20.7 mg of pantoprazole or the equivalent 23.4 mg of pantoprazole sodium standard. If you weighed 20.8 mg of pantoprazole standard, dissolve it in: (20.8 mg)/(0.383 mg/mL) = 54.3 mL of methanol. This makes the high standard solution concentration equal to 0.383 mg/mL. Low Standard:

The low limit is 85%; therefore the concentration of the low standard = (0.333 mg/mL X 0.85 = 0.283 mg/mL. Dilute 1.00 mL of high standard to 1.35 mL by adding 0.35 mL of methanol (1.15/0.85 = 1.35).

Spotting:

Spot on the 5 X 10 cm silica gel TLC aluminium plate with 3.00 μL aliquots as follows:

Left Spot low standard (85%) = $0.850 \mu g$

Center Spot 100% sample = $1.00 \mu g$

Right Spot high standard (115%) = $1.15 \mu g$

Development:


Mix 24.0 mL of ethyl acetate, 6.00 mL of methanol, and 12.0 mL of toluene. Develop the plate in a small glass chamber with approximately 20.0 mL of this solution until the solvent front reaches within 1 cm of the top of the TLC plate.

 $(R_f = 0.43)$

Detection:

UV:

Dry the plate and observe under ultraviolet light at 254 nm. Observe the intensities and the sizes of the spots.



Developed and tested by Kaitlin Nguyen and Joseph Sherma
Department of Chemistry, Lafayette College, Easton, PA, USA
October 2018