Sertraline HCl 50 mg Tablet

Structure:

Molecular Formula and Mass: C₁₇H₁₈Cl₃N – 342.688 **Category:** Selective Serotonin Reuptake Inhibitor

Sample:

Grind one tablet and dissolve in 15.0 mL of methanol. Shake for at least 10 min and filter. Final concentration of sample solutions = 3.33 mg/mL, which is the required concentration representing 100%.

Standards:

High Standard:

The high limit is 115%; therefore the concentration of the high standard = 3.33 mg/mL $\times 1.15 = 3.83$ mg/mL. Weigh approximately 38.3 mg of standard and dissolve it in 10.0 mL of methanol. If you weighed 38.2 mg of standard, dissolve it in: 38.2 mg $\div 3.83$ mg/mL = 9.97 mL of methanol. This makes the high standard solution concentration equal to 3.83 mg/mL, which is 115%.

Low Standard:

The low limit is 85%; therefore the concentration of the low standard = $3.33 \text{ mg/mL} \times 0.85 = 2.83 \text{ mg/mL}$. Dilute 1.70 mL of high standard to 2.30 mL by adding 0.60 mL of methanol. This gives a concentration of $3.83 \text{ mg/mL} \times 1.70 \text{ mL} \div 2.30 \text{ mL} = 2.83 \text{ mg/mL}$, which is 85.0%.

Spotting:

Spot on the 5×10 cm silica gel TLC aluminum plate with 3.00 µL aliquots as follows:

Left spot low standard (85%) = $8.49 \mu g$

Center Spot 100% sample = 9.99 µg

Right Spot high standard (115%) = 11.5 μ g

Development:


Mix 32.0 mL of toluene, 8.00 mL of ethyl acetate, 0.400 mL of ammonia, and 3.60 mL of ethanol. Develop the plate in a small glass chamber with approximately 20.0 mL of this solution until the solvent front reaches within 1 cm of the top of the TLC plate.

 $(R_f = 0.40)$

Detection:

UV:

Dry the plate and observe under ultraviolet light at 254 nm. Observe the intensities and the sizes of the spots.



Developed and tested by Bingsong Zeng and Joseph Sherma Department of Chemistry, Lafayette College, Easton, PA, USA April, 2018

Bingsong Zeng's EXCEL Scholar research was supported by Lafayette College EXCEL Scholars Program